

# hairotom

**Hairotom™** - a device for non-invasive collection of hair samples including follicular cells from mammals

## Abstract

In accordance with the 3 Rs principle (to replace, reduce and refine) animal models in biomedical research, a new approach is presented for sampling and analyzing hair follicles in various experimental settings. **Hairotom™** is a convenient device operable by a single person for the non-invasive collection of hair follicular cells from various mammals. The amount of collected biological material is sufficient to replace stressful and painful biopsies. Moreover, the main components of collected hair follicles are live cells of epithelial origin, which are highly relevant for studying ageing-related pathologies, including cancer. Obtained hair follicular cells are amenable for genotyping, quantitative PCR, and quantitative immunofluorescence using standard or slightly modified experimental protocols.


## Objectives

At present, invasive and non-invasive methods are used for collecting hair samples. The invasive procedure involves excising the tissue near the hair's root and processing the extracted tissue further. Invasive surgery always poses a risk of infection, permanent scarring and discomfort and with some hair types, for example, eyelashes and eyebrows, application of this method is very difficult.

Non-invasive methods involve extracting the hair in a way that the root surrounded by follicular cells is included. For this purpose, various types of ordinary tweezers, forceps and hemostats are used. The disadvantage of such regular metal instruments is the risk of breaking the hair in the area weakened by the grip. Another important factor in hair extraction is how easily the extraction device can grip many hairs. Finally, further manipulation with the extracted hairs introduces multiple practical scrutinies reflecting such samples' small size and general vulnerability.

## The technology

**Hairotom™** provides a new future in sampling and processing in translational research and routine applications, with a broad range of ethical and logistic advantages over currently used biopsy-based approaches. The device was optimized and validated for extraction of hairs from mammals to obtain complex hair samples allowing further analysis. The method and device is patent protected in EU and USA (WO 2014/086324, EP 2928382, US 14/649785)



The extraction device consists of a hand-held pistol connected to a vacuum source (e.g. standard vacuum cleaner). Another integral part of the extraction device is a double-shelled disposable forceps comprised of two hollow truncated telescopic cones. The cones are manufactured from a soft, inert and sterilization-compatible material resistant to the organic solvents.


## Method Benefits

- Non-invasive
- Easy to use
- Low costs
- Multiple hairs per sample
- Compatible with various mammalian hairs
- Disposable forceps
- Sterile conditions compatible

## Applications

- Genotyping
- Treatment response
- Biomarkers studies
- Basic research
- Biodosimetry
- Alopecia research
- Hair-associated microbiome analysis
- Derivation of primary cell lines


## Method overview




Device assembly is fast and easy, two plastic conical forceps are attached, and the whole pistol is connected to the vacuum source



A single person can perform the collection process from mice. The vacuum sucks the hairs between the forceps squeezed by the thumb operated piston.



Mice's hairs fixed in the forceps are easily pulled out. Each mouse hair includes a root area with approximately fifty follicular cells (microscopic inset image).





The firm grip of forceps also allows the withdrawal of deeply rooted hairs such as human hairs.




Collected samples can be processed in multiple ways, including direct fixation or further ex-vivo cultivations. The forceps are inert (fully compatible with organic fixatives), sterile, toxin- and pyrogen-free.

## Examples of Analysis



Collected hair samples can be used as a source of DNA and RNA. A) Results of genotyping of a transgenic mouse. The classical invasive approach of obtaining mouse DNA (tail clipping) is compared to the collected hairs' DNA. B) Assessment of mRNA quality isolated from hair follicular cells by miRNeasy mini kit (Qiagen) measured on Agilent RNA Pico Chip.



Various microscopic analyses were performed on the collected hair samples. A) 3D reconstruction of the sample using light-sheet microscopy (DAPI stained cell nuclei) B) Paraffin section of the hair stained by haematoxylin-eosin C) Immunofluorescence analysis of mouse follicular cells from gamma rays irradiated mouse stained for DNA (DAPI - blue) and marker of DNA damage (γ-H2AX - green). Hair body autofluorescence is visible in the red channel. (for more practical examples see also [www.aging-us.com/article/203744/text](http://www.aging-us.com/article/203744/text))



**AnLab, s.r.o.**  
Videňská 1083  
142 20 Prague 4  
Czech Republic

Phone: +420 261 711 667  
Mobile: +420 776 333 007  
Fax: +420 261 711 719  
E-mail: [info@anlab.cz](mailto:info@anlab.cz)

[www.anlab.eu](http://www.anlab.eu)  
[www.hairotom.com](http://www.hairotom.com)